Reference Sheet



This reference sheet provides an overview of all the keyboard shortcuts that are available and some help full tips for working with Cytosketch. You can also find the short cuts in the menu next to each item and in the tool tip text.

General Short Cuts

		-		
Ctrl +	N	曲	Create a new canvas	

Image Processing Short Cuts

- M Pages the images that are selected
- Splits an image into the separate channels
- Activates the area background tool
- Activates the area brightness tool
- Ctrl + T Auto Background
- Ctrl + B Auto Brightness
- Ctrl + G Auto Adjust
- Ctrl + R Resets selected images to original values

Short Cuts to Change the Zoom

- 700m in
- Zoom out
- Zooms image to 1 to 1
- Zoom to image
- Zoom to the selection
- Zoom to figure
- 2 Zoom 1 to 1

Short Cuts to Turn Channels On and Off

INUIT	nber	SHOV	VS C	me	cnanne	Office	
	_						

Shift + Number Toggles the channel on and off

To show all channels again

Channels

Number:

Channel:

Tips and Tricks

This section describes some useful tips and tricks for working with Cytosketch.

Tools

Keep shift pressed to use the tool multiple times. To exit the tool, press Esc or click the Select button \Im

Shift can be used with the following tools:

- Enumeration
- Pan
- Background
- Brightness

Un-proportionally Resizing

To un-proportionally resizing an image, keep Ctrl pressed while resizing the image.

Legend

Right click on the image to display the image legend.